

Mobile Response Unit

NON-EMS MOBILE MEDICAL VAN

Patent Pending | Not a Certified Ambulance



- ✓ Custom aluminum cabinetry layouts to best suit your organization's needs
- ✓ Designated storage spaces for special equipment such as stair chairs, spine boards, oxygen bottles, and more
- ✓ AutoFloor track system gives you full flexibility for stretcher and attendant seating layouts
- ✓ Clean, bright medical environment provides a comfortable staff and patient experience

A **mobile treatment space** and **transportation solution** for non-urgent medical needs.

RAM PROMASTER CHASSIS

Ram Promaster 159" Wheelbase Model
Front Wheel Drive
3.6L gasoline engine




MOBILE RESPONSE VAN CABINETRY

- Custom cabinetry layouts that work for your organization's medical service needs
- Painted aluminum finish provides a clean, bright medical environment
- Secure cabinet fronts give you easy access to medical supplies, medication, and equipment
- All surfaces are wipe-clean and easy to disinfect

STRETCHER AND ATTENDANT SEATING

- Layouts can include stretcher securement and/or attendant seating that secures right into the floor
- AutoFloor track system allows you to add and remove seats and change seat layouts on the go
- To move seats, simply lift a lever and wheel to desired location
- Attendant seats are grey wipe-clean transit grade leatherette



Easy access cabinets with durable, transparent slide and click fronts



Rear Heat & Air Conditioning controls in rear compartment. Air is directed at the discretion of staff for ultimate comfort. Control panel includes power outlets and lighting controls.



Contact us to discuss your mobile response van needs.

You can choose from a variety of standard layouts, or work with us to design a van that's fully customized for your needs.

MEDICAL EQUIPMENT STORAGE

All medical equipment has a designated storage space to make it easily accessible.
Additional options like a folding desk space or fridge are available.



First responder kit bag storage area.



Oxygen bottle holders, and stretcher securement.



Medical Glove Storage



Stretcher securement for major brands including Ferno and Stryker



Cab divider available with or without window.
Available with passthrough and sliding door options.



Serrated anti-slip running boards prevent ice and snow buildup



NON-EMERGENCY MOBILE MEDICAL VAN APPLICATIONS

AVAN Mobility's Mobile Response/Non-Emergency Mobile Medical Van is adaptable across many sectors and use cases, including:

- Native Reserves
- NEMT services, healthcare providers, outreach support organizations
- Housing societies
- Homeless and addictions programs

Contact us today to discuss your mobile medical van needs.